00657 Futura AD - Instrucao de Uso

190mm

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CE FUTURA AD Addition Silicone Impression Material English PROFESSIONAL USE ONLY

Composition:
FUTURA AD contains: Polyvinylsiloxane, Silicon Dioxide, Hydrocarbons, Methylhydrogenosiloxane, Lilac Coloring, Yellow Coloring, Blue Coloring and Organoplatinic Complex.

Indications:
FUTURA AD PUTTY base, lilac color, is indicated, after the appropriate mixture with FUTURA AD PUTTY catalyst, white color, as the first stage of molding. The indicated technique is the double impression and double mixture. It is also indicated for simultaneous molding with the Regular or Light Fluid.

FUTURA AD REGULAR FLUID base, yellow color, is indicated, after the appropriate mixture with FUTURA AD REGULAR FLUID catalyst, as a precision material on the second stage of molding. The indicated technique is the double impression and double mixture. It can also be used for simultaneous molding with the Putty. This consistency may also be used on the monophasic technique

FUTURA AD LIGHT FLUID base, blue color, is indicated, after the appropriate mixture with FUTURA AD LIGHT FLUID catalyst, as precision material on the second stage of molding. The indicated technique is the double impression and double mixture. It can also be used for simultaneous molding with the Putty.

FUTURA AD, provides:
Short Net Setting Time;
Easy handling, as it allows a self-mixture on the products with self-mixture cartridges.
Appropriate break endurance
High details precision
Inexpressive distortion after removal
Dimensionally steady even after one week.

Presentations:
•Kit containing 2 double cartridges with 75g(50 mL) each, of FUTURA AD REGULAR FLUID (Base + Catalyst) + 1 pot of 389g (250 mL) of FUTURA AD PUTTY Base + 1 pot of 380g (250 mL) of FUTURA AD PUTTY Catalyst + tips (mixers and intraorals) + 2 dosing spoons.
•Kit containing 2 double cartridges with 74g (50 mL) each, of FUTURA AD LIGHT FLUID (Base + Catalyst) + 1 pot of 389g (250 mL) of FUTURA AD PUTTY Base + 1 pot of 380g (250 mL) of FUTURA AD PUTTY Catalyst (mixers and intraorals) + 2 dosing spoons.
•Kit containing 1 double cartridge with 75g (50 mL) of FUTURA AD REGULAR FLUID (Base + Catalyst) + 1 double cartridge with 74g (50 mL) of FUTURA AD LIGHT FLUID (Base + Catalyst) + 1 pot of 389g(250 mL) of FUTURA AD PUTTY Base + 1 pot of 380g (250 mL) of FUTURA AD PUTTY Catalyst + tips (mixers and intraorals) + 2 dosing spoons.
•Kit containing 1 double cartridge with 74g (50 mL) of FUTURA AD LIGHT FLUID (Base + Catalyst) + 1 pot of 156g (100 mL) of FUTURA AD PUTTY Base + 1 pot of 152g (100 mL) of FUTURA AD PUTTY Catalyst + tips (mixers and intraorals), + 2 dosing spoons.
•1 pot of 389g (250 mL) of FUTURA AD PUTTY Base + 1 pot of 380g (250 mL) FUTURA AD PUTTY Catalyst + 2 dosing spoons.

•1 double cartridge with 74g (50 mL) of FUTURA AD LIGHT FLUID (Base + Catalyst) + tips (mixers and intraorals).
•1 double cartridge with 75g (50 mL) of FUTURA AD REGULAR FLUID (Base + Catalyst),+ tips (mixers and intraorals).
•12 mixer tips + 12 intraorals tips.

Storage and Conservation:
The molds with FUTURA AD should be preserved at room temperature. The FUTURA AD should be stored in a temperature between 2° C and 30° C. Keep the pots, the double cartridges and the tubes tightly closed after usage. For better conservation keep it in refrigerator. The product expiry date is 2 years after its manufacturing. Do not use this product if shelf life had expired.

Instructions for Use:
Working and Net Setting Time of FUTURA AD: PUTTY, REGULAR FLUID and LIGHT FLUID.

FUTURA AD	
Mixing Time*	30 seconds
Total Working Time*	2 minutes
Staying in oral cavity	3 min. and 30 sec.
Net Setting Time*	5 min. and 30 sec.

*High Temperatures accelerate the hardening and vice-versa.


Mold: Preparation and Adhesives:
Select the mold. All kinds of hard mold may be used, specially the stainless steel ones. It may also be used the triple mold with the FUTURA AD REGULAR FLUID, individuals acrylic molds or cut plastic molds for implant molding.
The adherence of the material to the mold represents an important phase to avoid mistakes, mostly in the phase when the mold is removed from the mouth. To assure a strong bond, it should be used only the Addition Silicon Adhesive.
Apply a fine layer of the adhesive on the mold and wait until it dries.

FUTURA AD PUTTY
As the Addition Silicon reacts when in contact with latex, do not use latex gloves when handling the FUTURA AD PUTTY.

Using the respective dosing spoons, retrieve not excessively and in equal portions, the necessary amount of FUTURA AD PUTTY base and PUTTY catalyst.
Mix identical proportions of Base and Catalyst, until you obtain a homogeneous colored mix, with no streaks (for 30 seconds).
Spread the formed paste over the mold. Before molding, make a relieve with PVC film or any other plastic, so that the elements to be molded plus a distal and mesial tooth save a space of 0.3 mm for FUTURA AD FLUID. The other elements should not be covered with this film.

Do not proceed with relieve in case you chose the simultaneous molding technique, which initiates with the handling of both consistencies at the same time. Await the complete Setting of the material and remove the mold from the mouth. Verify if the film kept the desired space. Otherwise, improve the relieve with the aid of a slightly cut milling machine.

FUTURA AD LIGHT OR REGULAR FLUID Cartridge Version:

- 1- Push up the small lever (A) behind the dispenser and pull the piston (B) totally from inside of the dispenser.
- 2- Raise the lock (C) over the dispenser, insert the cartridge. Align the crack on the cartridge base with the dispenser crack. Push hardly the lock to fasten the cartridge.
- 3- Before screwing the mixer tip, remove the cartridge cover (D).
- 4- Align the fitting of the mixer tip (E) and the casual intraoral tip (G) with the correspondent projection of the cartridge. Lock turning 90° clockwise.
- 5- Dispense the material through the mixer point pressuring the trigger (F) of the dispenser.

Attention: Before assembling the tip, ensure that the two components (base and catalyst) flow out evenly by applying a light pressure on the dispenser lever and extruding a small amount of material, which must be removed.

Tubes Version:
Dispense the necessary amount of FUTURA AD LIGHT or REGULAR FLUID base in length of the scale in a mixture block and then, an equal amount in length of FUTURA AD LIGHT or REGULAR FLUID catalyst.
Using the mixing block or glass plate, mixture energetically with an spatula eliminating the air bubbles contained in the mixture, until getting a homogenous color. The ideal mixing time is of about 30 seconds. Introduce the manipulated material in a molding syringe. After the syringe filling the rest of material is taken to the heavy base.

Use:
Double Mixture Technique:
After the gingival receding, if necessary, dry the elements to be molded. After removal of the retractor wire (if the two-wire technique is employed, the thinnest one may be kept on the gingival ridge), apply FUTURA AD LIGHT or REGULAR FLUID, first on the prepared elements and then on the other elements.
Locate the mold with FUTURA AD PUTTY immediately over the arch or hemi-arch to be molded. Wait for the complete cure of the material and remove the mold, checking if occurred any failure. Remove the retractor wire (for the two-wire technique) and proceed the next steps of the technique such as the agonist molding, the interocclusal registration.

Triple Molding Technique (for FUTURA AD REGULAR FLUID only)
After the gingival receding, if necessary, dry the elements to be molded. After removal of the retractor wire (if the two-wire technique is employed, the thinnest one may be kept on the gingival ridge), apply REGULAR FLUID, first on the prepared elements and then on the other elements.

Simultaneously, dispense the same material over the triple mold previously conditioned with adhesive for addition silicon.
Take the mold into the mouth and request the patient to apply a light pressure. Wait for the complete cure of the material and remove the mold, checking if occurred any failure.
Remove the retractor wire (for the two-wire technique)

Cleaning and Disinfection of the Molding:
After having been well rinsed under running water, over 1 minute, the mold can be disinfected immediately.
Casting Models:
ATTENTION: Wait for at least 1 minute after the molding to proceed with the casting models.
The casting models can be performed immediately after the disinfection and up to two weeks later, provided that the molding is kept at room temperature and the time of at least 1 hour after the molding is respected.
For a best quality of the mold, cast model type IV.
The use of addition silicon allows a second casting model.

Precautions:
Avoid the use of FUTURA AD in patients with allergy history to any of the components of the formula. In the case of allergic reactions, suspend the use immediately and tell the patient to seek medical advice.
Avoid the contact with the skin and the eyes with the catalyst. In the case of contact with the eyes, wash thoroughly with water and seek medical device. In case of contact with the skin, wash with water and soap. In case of accidental ingestion, seek medical assistance immediately. Avoid direct contact among the product and latex gloves, retraction wires or other preparations, once these can inhibit the product setting. It is being advised that the hands be very well washed to eliminate any pollutant trace.
The mentioned proportions should be respected to obtain the molds with correct setting. A bigger dose of the catalyst does not accelerate the time of net setting.
Close the packages immediately after its use and do not change the covers or the dosing spoons.
The polymethylsiloxanes are chemically resistant therefore, avoid spotting fabrics.
Leftovers or residues of used materials should not be reutilized.
They must be disposed as the local law in force.

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