
00663 Clonage - Instrucao de Uso PEI

Instructions:

	CLONAGE PUTTY	CLONAGE FLUID
Mixing time	30 sec	20 sec
Total Working Time*	1 min 15 sec	1 min 30 sec
Time in oral cavity	2 min 15 sec	2 min
Hardening time (fast set)*	3 min 30 sec	3 min 30 sec
Strain in compression (min - max)	2% - 5%	2% - 8%
Recovery from deformation	> 98%	> 99%
Dimensional stability	< -0,2%	< -0,7%

The times are counted from the beginning of the mixing at 23°C.

High temperatures accelerate the hardening and vice-versa. Excessive catalyst dosages do not accelerate the hardening time.

1st Step

Select a rigid mold, stainless steel ones are preferred.

Using the measuring scoop, withdraw (without excesses) the necessary amount of **CLONAGE PUTTY** and spread it out over your hand. Mark the **CLONAGE PUTTY** on the hand using the measuring scoop. For each dose of the material, apply a straight strip of **CLONAGE CATALYST** in the length and another across the material (forming a cross inside the edges of the mark). See Figure 1.

Fold the paste over itself and manipulate energetically for about 30 seconds, until the paste becomes homogenous.

Spread the paste over the mold. Before molding, make a relief with a PVC film or any other plastic on the more distal and mesial teeth to be molded, so that a space of about 0.3 mm for **CLONAGE FLUID** is left. The other elements should not be covered by this film.

Seat loaded tray. Wait until the complete setting of the material and remove the impression from the mouth. Check if the film kept the desired space. If not, improve the relief with a fine cut bur.

2nd Step

After gingival receding, if necessary, dry the elements to be molded.

Dispense **CLONAGE FLUID** in a mixing block or plate and an equal amount of **CLONAGE CATALYST**, with the same length of the **CLONAGE FLUID**. The thickness must be proportional to the diameter of the orifice of each tube. See Fig. 2.

Mix energetically with a spatula eliminating air bubbles entrapped in the mixture, until a homogenous pink color is obtained.

The ideal mixing time is about 20 seconds. Fill the impression syringe with the manipulated material and dispense the rest of the material over the putty impression.

Remove the retractor cord and inject **CLONAGE FLUID** with the syringe, first on the prepared elements and then on the other elements.

Seat the loaded tray. Wait until the complete setting of the material and remove the impression from the mouth.

Casting

Before pouring gypsum, wait 20 minutes for the elastic recovery of the impression. During this time, the impression should be disinfected, rinsed with running water for 1 minute and dried.

The impression should be poured from 20 minutes to 1 hour after hardening. For a high model quality, pour type IV gypsum, wait from 50 minutes to 1 hour and remove the impression from the gypsum.

Precautions:

Avoid the use of **CLONAGE** in patients with history of allergy to any of the components in the formula. In the case of allergy, suspend the use and seek medical advice.

Avoid the contact of the skin and the eyes with the catalyst. In the case of contact with the eyes, wash thoroughly with water and seek medical advice.

In case of contact with the skin, wash with water and soap.

In the case of accidental ingestion, seek medical advice immediately.

This product must be used only by dentists and for the intended use described in the Instructions for Use.

Product wastes must not be reused. They must be disposed as the local law in force.

Fig.1

Fig.2

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CLONAGE Condensation Silicone Dental Impression Material English

CLONAGE is a high quality and fast set condensation silicone developed to take dental impressions that require high precision and definition of details.

CLONAGE PUTTY presents:

- Unique rheological property, where the fluidity varies in function of exerted effort;
- Elimination of patient nausea and discomfort;
- High softness during application, combined with better rigidity after setting time.

CLONAGE FLUID has a unique formulation, developed to maximize the dimensional accuracy of the impression. CLONAGE FLUID presents:

- Better penetration in subgingival areas, assuring better reproduction of details;
- Hydrophilicity, enabling impressions in the presence of humidity;
- Homogeneity, forming a uniform layer over CLONAGE PUTTY;
- Excellent color contrast with CLONAGE PUTTY;
- Easy removal from mouth, after setting.

Composition:

CLONAGE PUTTY	contains: Polydimethylsiloxane.
CLONAGE FLUID	contains: Polydimethylsiloxane.
CLONAGE CATALYST	contains: Alkyl Tin Siloxane.

Applications:

After adequate mixture with CLONAGE CATALYST, CLONAGE PUTTY is recommended to take the first impression in the two-step technique (tray material) and CLONAGE FLUID is recommended to take the precise impression in the two step technique (wash material).

Presentations:

Box containing 1 pot with 1000 g of CLONAGE PUTTY.

Box containing 1 tube with 120 g of CLONAGE FLUID.

Box containing 1 tube with 50 g of CLONAGE CATALYST.

Box containing 1 pot with 500 g of CLONAGE PUTTY.

Box containing 1 tube with 60 g of CLONAGE FLUID.

Box containing 1 tube with 25 g of CLONAGE CATALYST.

Storage and Conservation:

CLONAGE must be kept away from excessive heat (avoid temperatures higher than 40°C). Keep the pot and the tubes tightly closed after use. The product has a shelf life of 2 years. Do not use this product if shelf life had expired.